

DENTRIX

Dentrix Learning Edition

SYSTEM REQUIREMENTS

PUBLICATION DATE

August 2020

COPYRIGHT

© 1987-2020 Henry Schein One. All rights reserved. No part of this publication may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any language in any form by any means without the prior written permission of Henry Schein One.

SOFTWARE LICENSE NOTICE

Your license agreement with Henry Schein One, which is included with the product, specifies the permitted and prohibited uses of the product. Any unauthorized duplication or use of Dentrix by Henry Schein One in whole or in part, in print, or in any other storage and retrieval system is forbidden.

LICENSES AND TRADEMARKS

Henry Schein One, the Henry Schein One logo, and Dentrix are registered trademarks of Henry Schein One. Microsoft and Windows are registered trademarks of Microsoft Corporation.

Dentrix Learning Edition System Requirements

The System Requirements describe minimum standards for using Dentrix Learning Edition. Requirements can vary significantly depending on workload and other software being used. Exceeding the minimum standards may result in better system performance.

Over time, system requirements change. For the latest system requirements, visit learningedition.dentrix.com.

| | |
|----------------------------------|-------------------------------|
| Operating System | Windows* 7 Sp1, 8.1, 10 |
| Memory | 8 GB RAM |
| CPU | 4 cores at 2.4 GHz |
| Local Drive Install Space | 40 GB total, 5 GB on C: drive |
| Monitor | 1280x1024 |

*Dentrix Learning Edition supports any edition of the Windows workstation operating system such as Home, Pro, Enterprise, and so forth.

Notes & Additional Recommendations

TERMINAL SERVICES

Thin client setups, such as Terminal Services and Citrix, are currently not supported and should not be used with Dentrix.

HARD DRIVE

The disk space needed for Dentrix depends on the size of the practice and the amount of data that will be stored. For best performance, we recommend more than 10 percent free space on all physical drives.

NETWORKS

High-speed internet connectivity is recommended for access to software updates and all available services.

PRINTERS

Choose a printer based on your practice needs. Henry Schein One cannot guarantee that all printers will be completely compatible with Dentrix. We recommend that you use PCL5 printer drivers with all printing equipment.

MICROSOFT OFFICE

Printing form letters and labels from Dentrix requires a current Microsoft-supported version of Microsoft Word with letter merge functionality to be installed.

ANTIVIRUS SOFTWARE

Antivirus software is recommended. When configuring antivirus software, consider the inherent risks in excluding certain items from the protection afforded by antivirus software.